

Installation Guide

Adaptor Plate for GE Corometrics 170 Series
(Refer to qualified personnel)

Purpose:

To describe the installation of the Adaptor Plate onto the GE Corometrics 170 Series device and placement onto GCX Fetal Monitoring Cart



Parts Reference

The following parts and hardware are included in this installation kit:

Item #	Description	Qty
1	Adaptor Plate	1
2	4-40 x 1/2" Flat Head Machine Screw (FHMS)	4

Tools Required:

- Phillips Screwdriver




Removal of the Rubber Feet on the Device

1. Remove the four feet from the bottom of the device.



Installation of the Adaptor Plate onto the GE Corometrics 170 Series Device

1. Fasten the Adaptor Plate with four (4) 4-40 x 1/2" Flat Head Machine Screws (FHMS) onto the GE Corometrics 170 Series Device with a Phillips screwdriver as shown below.



Installation of the Adaptor Plate onto Fetal Monitoring Cart Top Surface

1. Lower the Adaptor Plate onto the Tray and secure it by placing the four Locating Pin on the Tray inside the four Locating Holes on the Adaptor Plate.



Routine Maintenance

Periodically inspect all fasteners associated with the mounting assembly. Tighten or Adjust as necessary for optimal operation and safety.

Cleaning the Mounting Assembly

1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone or trichloroethylene.
3. Steel wool or other abrasive material should never be used.
4. Damaged caused by the use of unapproved substances or processes will not be warranted. We recommend testing of any cleaning solution on a small area of the mounting assembly that is not visible to verify compatibility.
5. Never submerge or allow liquids to enter the mounting assembly. Wipe any cleaning agents off the mounting assembly immediately using a water-dampened cloth. Dry the assembly thoroughly after cleaning.

CAUTION: GCX makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted devices or accessory equipment, refer to the specific instructions delivered with those products.