

## GCX Mounting Assembly Installation Guide

### Low Profile Countertop/Horizontal Channel Mount

*The purpose of this guide is to:*


1. Describe the procedure for mounting the mounting assembly in a horizontal channel.
2. Describe the procedure for mounting the mounting assembly on a counter top.
3. Describe procedures for swiveling and tilting the monitor and adjusting tilt tension.

### Parts Reference

The following parts and hardware are included for this installation procedure (hardware not shown):

Item #	Description	Qty
1	Low Profile Countertop/Horizontal Channel Mount	1
*2	1/4-20 x 1-3/4" Carriage Bolt	4
*3	1/4-20 x 3/4" Carriage Bolt	4
*4	1-1/4" (OD) Flat Washer	4
*5	1/4-20 Kep Nut	4
6	5/32" Hex Key	1

*\*Required if mounting on a counter top.*



### Tools Required:


- Drill with 17/64" drill bit (not provided)
- 7/16" wrench (not provided)
- 5/32" hex key (provided)
- Phillips Screwdriver (not provided)



**WARNING:** It shall be the responsibility of the installer to determine whether the counter top surface will support the weight of the Mount and any mounted instruments or accessories.

### Installing the Mounting Assembly in a Horizontal Channel


1. Slide the Base Plate into the channel and move to desired mounting position.
2. Tighten two (2) thumb screws to secure position in channel.



### Installing the Mounting Assembly on a Counter Top


**Installation Note:** Before installing the mounting assembly on a countertop, make sure that thumb screw threads are not protruding from the bottom of mounting assembly.

1. Using the Base Plate as a template, drill four (4) holes in the counter top with a 17/64" drill bit.
2. Fasten Base Plate to counter top using four (4) 1/4-20 carriage bolts, four (4) 1-1/4" flat washers, and four (4) 1/4-20 kep nuts. Tighten nuts using a 7/16" wrench (not provided).




## Adjusting Swivel and Tilt Tension

Adjust **tilt tension** by tightening or loosening the Tilt Tension screw on the top of the Mount. A 5/32" hex wrench is provided for this adjustment.



Adjust **swivel tension** by tightening or loosening the swivel tension screw on the side of the Mount. A Phillips screwdriver (not provided) is used for this adjustment.



## Routine Maintenance

Periodically inspect all fasteners associated with the mounting system. Tighten as necessary for optimal operation and safety.

## Cleaning the Mounting System

1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone and trichloroethylene.
3. Do not use steel wool or other abrasive material to clean the mounting assembly.
4. Damage caused by the use of unapproved substances or processes will not be covered by warranty. We recommend testing of any cleaning solution on a small area of the mounting assembly that is not visible to verify compatibility.
5. Never submerge or allow liquids to enter the mounting assembly. Wipe any cleaning agents off of the mounting assembly immediately using a water-dampened cloth. Dry the mounting assembly thoroughly after cleaning.

**CAUTION:** GCX makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted instruments or accessory equipment, refer to the specific instructions delivered with those products.