Mindray V Series VHM Wall Mount

Installation & Operation Manual

The purpose of this manual is to describe general installation, operation, and adjustment procedures for the Mindray V Series VHM Wall Mount. Please read this manual before installing or using this product.

Table of Contents

1.0. Installation, Operation, and Maintenance Warnings.....	2
2.0 Installing the VHM Arm in the Channel.....	3
3.0 Mounting the Instrument on the VHM Arm.....	4
4.0 Operating and Adjusting the VHM Arm.....	4
5.0 Cable Management.....	9
6.0 Routine Maintenance.....	9
7.0 Cleaning the Mounting Assembly.....	9
8.0 Trouble shooting the VHM arm	10

1.0. Installation, Operation, and Maintenance Warnings

This section contains warnings regarding the installation, operation, and maintenance of the VHM Arm. This section must be read in its entirety before installing, operating, and maintaining the VHM Arm. Failure to follow these warnings may result in damage to equipment or injury to personnel.

Warnings

- Do not position VHM Arm or mounted instrument above a patient. Note that the VHM ARM has a wide range of motion both up/down and side to side. Please consider carefully the instrument being mounted and the proximity of the mounting assembly to other equipment, hospital personnel, and the patient. GCX recommends that the hospital's risk management personnel verify that the application is appropriate prior to installation and use of the VHM Arm.
- If VHM Arm is wall mounted, verify that the channel has been installed and approved in accordance with channel installation guide.
- *Ensure that the weight of the instrument being mounted does not exceed the load rating of the VHM Arm. Check the bottom surface of the Arm for "Maximum Load" label (see 2.1 for label location). Mounted instrument must be within "Optimal Weight Range" for best performance.*
- The Height Locking Lever or Knob must be engaged when removing a mounted instrument from the VHM Arm.
- Stand to the side of the VHM Arm and mounted instrument and use caution when disengaging the Height Locking Knob. The total mounted load may have changed due to the addition or removal of equipment (other than primary instrument) while the Height Locking Knob was engaged. A change in the total load can cause a sudden downward or upward movement of the Arm when the Height Locking Knob is disengaged.
- The mounted device may move suddenly due to normal wear or improper adjustment of the tilt and swivel functions (see Sections 4.0 and 4.1) or ultimately, gas spring end of life. The VHM Arm must be inspected and adjusted at least two times per year. This inspection must include the following steps:
 - 1) With the device(s) mounted, release the lock and move the arm through its entire vertical range of motion. The load should maintain its position at every point in the travel of arm. If necessary, the counterbalance mechanism may be adjusted in accordance with section 4.0.
 - 2) Grasp the mounted device and swivel it from side to side. The device should swivel with some tension or resistance, not loosely. If necessary, the swivel tension may be adjusted in accordance with section 4.5) Grasp the Arm and pivot it from side to side at the channel attachment. The arm should pivot with some tension or resistance, not loosely. If necessary, the pivot tension may be adjusted in accordance with section 4.4.
 - 4) Grasp the mounted device and tilt it forward and back, through its entire range of motion. There should be enough tension or resistance in the tilt mechanism to prevent the device from tilting forward unexpectedly when in use. If necessary, the tilt tension may be adjusted in accordance with section 4.6.
 - 5) Inspect all remaining fasteners for looseness. Tighten as required for optimal operation and safety.
 - 6) The VHM Arm housing must never be disassembled by non-GCX personnel.

Failure to periodically inspect and adjust the Arm as instructed may result in damage to equipment or injury. If properly maintained, the VHM Arm should function properly for 5 to 10 years. If you are unable to adjust the arm or require service, contact GCX technical support at 800-228-2555.

Install your VHM Arm in the sequence presented in this manual (Sections 2.0, 3.0, and 4.0).


2.0 Installing the VHM Wall Mount in the Channel

2.1 VHM Arm with Spring Lock

This type of VHM Arm is designed for use with transport monitors only. The Locking Knob is installed at the factory. Ensure Locking Knob is in locked position before attempting to mount Arm in channel.


Installation Note: For proper function of the locking mechanism, the arm must be properly adjusted to counterbalance the weight of the mounted instrument (see section 4.1).

Note: Locking Knob operates by pulling and releasing Knob.


2.2 While supporting the bottom of the Arm, guide the Slide into the top of the channel (below left) and move Arm to required mounting position.

2.3 Using the 1/8" hex wrench provided, tighten four (4) set screws in Slide to secure position of Arm (below right).

3.0 Mounting the Docking Station on the VHM Wall Mount

3.1 Pull the plunger at the front of the mounting plate. Slide the Mindray V Dock mounting adapter into the mounting plate until the plunger snaps into the clearance hole in the adapter. Tighten the (4) Nylon Screws on the underside of the plate.

4.0 Operating and Adjusting the VHM Wall Mount

Check pivot, tilt, and front swivel mechanisms for proper tension. Refer to the applicable section and follow the adjustment procedure.

4.1 Adjusting Counterbalance

WARNING: Do not attempt counterbalance adjustment unless the instrument and accessories are mounted on the Arm. Use caution while performing this procedure.


The VHM Arm must be adjusted to counterbalance the mounted instrument throughout the Arm's vertical range of motion. When properly counterbalanced, the VHM Arm will maintain its height without engagement of the Height Locking Knob. Because instrument weights vary, adjustment will typically be required to achieve optimal performance. A 1/2" [13mm] socket wrench is provided for this critical adjustment.

To adjust counterbalance:

4.1.1. Grasp the Arm behind the mounted instrument(s) and carefully move the Arm to a horizontal position that allows access to the Adjuster Cover (below left). Lock Arm in horizontal position.

4.1.2. Open the Adjuster Cover by inserting a flat blade screwdriver in the slot at the rear of the Cover and prying upward. If necessary unlock and readjust the Arm until the Counterbalance Adjuster bolt is accessible through the Adjuster Cover (below center). Lock Arm when Counterbalance Adjuster is accessible.

4.1.3. With the 1/2" [13mm] socket wrench on the Adjuster, carefully unlock the arm. Turn the Adjuster *counterclockwise* (CCW) to increase counterbalance force, or *clockwise* (CW) to decrease counterbalance force. Counterbalance is correctly adjusted when the mounted instrument can be moved up or down with minimal force and does not rise or fall after releasing the Arm.

4.2 Adjusting Height of the Mounted Instrument

WARNING


- Do not attempt to adjust height without device(s) mounted on Arm.

- Stand to the side of the VHM Arm and mounted instrument and use caution when disengaging the Height Locking Knob. The total mounted load may have changed due to the addition or removal of equipment (other than primary instrument) while the Height Locking Knob was engaged. A change in the total load can cause a sudden downward or upward movement of the Arm when the Height Locking Knob is disengaged.

4.2.1. Grasp mounted instrument and move to desired height. When Arm is counterbalanced correctly (sec 4.1), it will maintain position when adjusted to a new height.

Note: Disengagement/ re-engagement of the Height Locking Knob is necessary during vertical adjustment of mounted instrument (see **WARNING** above).

 Use Caution when disengaging the Height Locking Knob. The mounted load may have increased (addition of equipment) or decreased (removal of equipment) while the Height Locking Knob was engaged, causing sudden downward or upward movement of arm.


 **Instrument Must Be Mounted**

4.3 Positioning VHM Arm in Channel

Caution: Do not attempt to remove the Arm from the Channel while Arm is loaded with an instrument.


- 4.3.1. Lock the Arm in highest vertical position using Height Locking Knob.
- 4.3.2. Using the 1/8" hex wrench provided, loosen four (4) set screws in Slide.
- 4.3.3. Relieve the weight of the mounted device(s) by lifting against the underside of the Arm, near the instrument. Simultaneously, push up or pull down the opposite end of the Arm immediately in front of the Slide and move to desired height.
- 4.3.4. Tighten four (4) Set Screws to secure position in Channel.

4.4 Pivoting the VHM Arm and Adjusting Pivot Tension

4.4.1. To pivot Arm at Channel, simply push on side of Arm.

4.4.2. To adjust pivot tension, remove upper and lower bolt caps from Pivot Bolt. Using the 1/2" [13mm] socket wrench provided and a 1/2" [13mm] box or open-end wrench, tighten or loosen Pivot Bolt and Hex Nut to desired tension. Reinstall bolt caps on each end of Pivot Bolt.


4.5 Swiveling the Mounted Instrument and Adjusting Swivel Tension

4.5.1. To swivel mounted instrument, simply push or pull instrument.

4.5.2. To adjust swivel tension, tighten or loosen the Swivel Tension Nut using the 1/2" [13 mm] socket wrench provided.

Installation Note: It may be necessary to remove an existing down post from the swivel cup to gain access to the Swivel Tension Nut.

Warning: Swivel tension nut must be torqued to a minimum of 20 in-lbs [2.3 N-m].


4.6 Adjusting Tilt Tension and Tilting the Mounted Device

In most configurations, the mounted instrument may be tilted and tilt tension may be adjusted.

4.6.1. Adjust tilt tension by equally tightening or loosening two (2) Tension Adjustment Screws with provided 5/32" Hex Wrench.


4.6.2. To tilt mounted instrument, simply push or pull instrument.

Warning: Tilt tension screws must be torqued to a minimum of 35 in-lbs [4.0 N-m].

5.0 Cable Management

Two (2) Cable Guides are provided to facilitate routing of cables along the underside of the Arm. Squeeze the Cable Guide until its edges snap into the grooves in the bottom surface of the Arm. Leave some extra cable loose at the front and rear of the arm to prevent cable binding, connector damage, or Cable Guide damage.

6.0 Routine Maintenance

Periodically inspect all tilting, swiveling, pivoting, and mounting hardware. Tighten or adjust as necessary for optimal operation and safety.

7.0 Cleaning the Mounting Assembly

- 7.1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
- 7.2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone and trichloroethylene.
- 7.3. Steel wool or other abrasive material should *never* be used.
- 7.4. Damage caused by the use of unapproved substances or processes will not be warranted. We recommend testing any cleaning solution on a small area of the arm that is not visible, to verify compatibility.
- 7.5. Never submerge or allow liquids to enter the arm. Wipe any cleaning agents off of the arm immediately using a water-dampened cloth. Dry the arm thoroughly after cleaning.

CAUTION: GCX makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted instruments or accessory equipment, refer to the specific instructions delivered with those products.

8.0 Troubleshooting the VHM Arm

Symptom	Possible Cause	Remedy
Mounted instrument does not appear level or parallel to the floor.	Channel not plumb. Check with level.	Adjust Channel to plumb, or reinstall Channel.
	Weight of instrument not compatible with Load Rating of the Arm.	Mount instrument on arm with compatible Load Rating.
	Swivel hardware loose.	Adjust Swivel Nut (section 4.5).
	Pivot hardware loose.	Adjust Pivot Bolt (section 4.4).
	Mounting surface (e.g. wall, side of anesthesia machine, etc.) not structurally sound (does not hold mounting hardware).	Mounting surface must be reinforced or Channel must be relocated.
	Channel loose at mounting surface.	Check for plumb and tighten, or relocate (reinstall) Channel.
Instrument drifts up or down when the Arm is unlocked.	Arm not counterbalanced correctly for weight of the instrument.	Perform counterbalance adjustment per section 4.1.
	Weight of mounted instrument (load) not compatible with Load Rating of Arm.	Use arm with compatible Load Rating, and perform counterbalance adjustment per section 4.1.
Arm pivots too freely.	Pivot Bolt too loose.	Adjust Pivot Bolt (section 4.4).
Arm does not pivot easily.	Pivot Bolt too tight.	
Instrument swivels too freely.	Swivel Nut too loose.	Adjust Swivel Nut (section 4.5).
Instrument difficult to swivel.	Swivel Nut too tight.	
Instrument difficult to tilt.	Tilt Adjustment Lever too tight.	Adjust tilt and tilt tension per section 4.6.
	Excessive tilt tension.	
Instrument will not maintain tilt position.	Tilt Adjustment Lever too loose.	
	Insufficient tilt tension.	
Arm difficult to move up or down when unlocked.	Arm not counterbalanced correctly for weight of mounted instrument.	Perform counterbalance adjustment per section 4.1.
Arm inadvertently slides down Channel.	Set Screws (4) in Slide are loose.	Reposition Arm and tighten set screws.
Locking Knob (Arms with Spring Lock) will not pull out (release).	Arm not counterbalanced correctly for weight of mounted instrument.	Perform counterbalance adjustment per section 4.1.
	Weight of mounted instrument (load) not compatible with Load Rating of Arm.	
Counterbalance Bolt difficult to adjust.	Arm not locked in a horizontal position.	Reposition Arm and lock in horizontal position.

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