


Counter Top Mounting Solutions for Datex-Ohmeda Light Monitor



Light monitor mounted on a counter top mount, P/N: 8000726. Requires Datex-Ohmeda supplied mounting adapter, P/N: 890915.



Counter Top Mount

Datex-Ohmeda P/N: 8000726

Features:

- 5" (12.7 cm) height
- Full swivel capability
- Adjustable tilt for optimal viewing angle
- Durable anodized aluminum components

GCX® offers a complete range of mounting solutions. Please contact us or visit our web site for additional information.

GCX
CORPORATION

3875 Cypress Drive
Petaluma, CA 94954-5635
707.773.1100
800.228.2555 US/CAN
Fax 707.773.1180
www.gcx.com

GCX is an ISO9001 registered company.

8/00-2



GCX is a registered trademark of GCX Corporation.

© 2000