

# Patient Engagement Table

Meet patient engagement objectives with a versatile solution that integrates technology and involves patients in their own care.



## GENERAL FEATURES

- Low-profile, U-shaped base for bed clearance and ease of movement
- Seamless, durable synthetic laminate surface with recessed area to contain spills
- Height adjustable via convenient, easy-to-use hand controls
- Anodized aluminum construction for durability and ease of cleaning
- Integrated, industry-standard GCX channel profile to facilitate IT hardware integration
- Accessory track to integrate accessories (storage, USB charging module, medical-grade power strip)
- Handle for easy positioning with cradle for power cord
- VHM-T Variable Height Arm for tablet device with full motion adjustment for easy positioning
- Shown with optional tablet enclosure\*

## ACCESSORIES

### USB Charging Modules

Compatible with most tablets and mobile phones module with coil-type or straight hospital-grade cord set

### Mountable Tablet Enclosures

Medical-grade mountable tablet enclosures for Apple, Samsung and Microsoft tablets - tablets sold separately.

### Medical Grade Power Strip SR-0001-07 + RST-0007-26

Medical-grade four hospital-grade outlet power strip - 110V Only - UL1363A - Approved - for OB T column

### Utility Bin WM-0001-29

For patient or clinician storage; easily removed for cleaning

### Flip-Down Worksurface OBT-0002-27


For OB T Rear Column Accessory Track -- 14 x 8" / 35.6 x 20.3 cm

## TECHNICAL SPECIFICATIONS

- Tabletop material: Laminate with Protective Edge\*\*
- Tabletop dimensions, overall w/d: 19.25 x 32" / 48.9 x 81.3 cm
- Tabletop dimensions, recessed area w/d: 18 x 30" / 45.7 x 76.2 cm
- Height adjustment range: 14.75" / 37.5 cm
- Tabletop height, minimum to maximum: 27.75 to 42.5" / 70.5 to 108 cm
- Main column height: 46" / 116.8 cm
- Base clearance height: 2" / 5.1 cm
- Caster diameter: 1.4" / 3.6 cm

Five (5) year warranty  
Patent Pending

\* Tablet enclosure to secure and protect the tablet device - Sold separately  
\*\* Tabletop options available



GCX CHANNEL PROFILE



UP / DOWN  
RELEASE



HANDLE



MOUNTABLE TABLET  
ENCLOSURE  
(OPTIONAL)



USB CHARGING  
MODULE  
(OPTIONAL)

## TABLETOP COLOR OPTIONS



GLACIER WHITE



CLOUD GRAY



LIGHT MAPLE



VALENCIA



MAHOGANY