



**Device/Mount Compatibility List for GCX Part Number:**

**FLP-0008-38 Roll Stand Post Mount for 75mm VESA Data Capture Devices with Adjustable Tilt and Height**

**This list is current as of July 2015. All configurations have passed tip testing based on a data capture device and accessories weighing a maximum of 6.6 lbs/3kg. Contact GCX to obtain the latest information on a specific configuration not listed or for information if using a device greater in weight than 6.6 lbs/3kg.**

**Mindray/Datascope**

Accutor V on GCX Roll Stand Part Number: RS-0006-64

Accutor Plus on GCX Roll Stand Part Number: RS-0006-64

**Philips**

VS2 on GCX Roll Stand Part Number: PH-0058-62

VS2\*/VS2+/VS3/V4, VS<sub>i</sub>, VM4/6/8, VSV on GCX Roll Stand Part Number: PH-0058-60

**WelchAllyn**

VS Spot LXI on GCX Roll Stand Part Number: RS-0025-02

VS Spot LXI on WelchAllyn Spot Vital Signs Monitor Roll Stand

**Criticare**

506DN eQuality on GCX Roll Stand Part Number: CT-0027-60

**GE**

DASH 3000/4000 on GCX Roll Stand Part Number: RS-0006-64 or ME-0030-63