


## Information System Component Mobile Mounting Solutions – VHRS™



### VHRS Variable Height Roll Stand with Foot Pedal for Flat Panel Display (Max. 23 lbs./10.4 kg)

#### FLP-0001-68

- VESA® FDMI™ 75/100 mm compatible
- Weight range: 5 - 23 lbs./2.3 - 10.4 kg
- Minimum height: 34.5"/87 cm
- Maximum height: 45.5"/115.5 cm
- 11"/28 cm height adjustment range
- 21"/53.5 cm diameter aluminum base built-in 10 lbs./4.5 kg counterweight
- 4"/10.2 cm diameter dual wheel nylon casters (2 locking and 3 non-locking) for easy transition over uneven surfaces
- Foot pedal control for height adjustability
- Gas spring-assisted for easy height adjustment
- A wide variety of accessories are available for practical convenience and cable management

### VHRS Variable Height (61.5"/156.2 cm) Roll Stand with Foot Pedal for Flat Panel Display (23 lbs/10.4 kg Max.)

#### FLP-0001-97

- VESA® FDMI™ 75/100 mm compatible
- Weight range: 5 - 23 lbs/2.3 - 10.4 kg
- Minimum height: 50"/127 cm
- Maximum height: 61.5"/156.2 cm
- 11"/28 cm height adjustment range
- 25"/63.5 cm diameter aluminum base built-in 21 lbs/9.5 kg counterweight
- 4"/10.2 cm diameter dual wheel nylon casters (2 locking and 3 non-locking) for easy transition over uneven surfaces
- Foot pedal control for height adjustability
- Gas spring-assisted for easy height adjustment
- A wide variety of accessories are available for additional components and cable management



GCX is an ISO 9001 registered company (TUV-Essen)  
GCX, VHRC, VHRS, VHM, PRC, M-Series, PolyMount and PolyQuip are trademarks of GCX Corporation.

© Copyright 2010

MKT-0001-19 2/2010