

VHRC™ for Philips IntelliVue MX600/MX700/MX800 with Laptop Enclosure

VHRC Variable Height Roll Cart

RST-0007-04

- The integrated VHRC channel accepts all GCX channel mounting components
- Cart overall height 64.5"/163.8 cm
- 40"/101.6 cm vertical adjusting channel offers 16"/40.6 cm height adjustment
- Integrated/adjustable cable management
- 18" x 18"/45.7 cm x 45.7 cm base
- 5"/12.7 cm diameter casters (4 locking)
- 50 - 60 lbs/22.7 - 27.2 kg Weight Range

M-Series Flush Mount for Philips IntelliVue MX600/MX700/MX800

AG-0018-73

- Tilt and swivel adjustment
- Tilt mechanism maximum load: 30 lbs/13.6 kg
- Product weight: 1.5 lbs/.7 kg

Laptop Enclosure

RST-0007-30

- Allows the notebook to be secure and used with an external keyboard and mouse
- Vented enclosure with rear openings for cables
- Capacity for 17"/43.2 cm wide notebook
- Left/right accessory compartments
- 18.5"/47 cm x 10.5"/26.7 cm work surface forward of the display
- Front handle for ease of movement
- Max weight for enclosure: 12.5 lbs/5.7 kg
- Product weight: 13.2 lbs/6 kg

Slide-Out Keyboard

RST-0007-34

- Slide out keyboard tray with padded wrist rest
- Keyboard tray 20.8"/52.8 cm x 8"/20.3 cm with left/right slide out mouse trays and Mouse pads (2) with adhesive backing
- Product weight: 18.2 lbs/8.3 kg

Flush Mount for Philips MP2/X2, FMS or Power Supply

HP-0053-07

- Ideal for applications with limited space
- Interface adapter plate and fixed flush channel slide for use in a GCX channel
- Adjustable stop offers desired height positioning
- Power Supply and MP2/X2 picture not shown

Note: Docking station for MP2/X2 is not included and must be ordered from Philips.

FMS Utility Hooks

HP-0053-08

- Utility Hooks For Philips Satellite Rack/FMS

Counterweight Kit

RST-0007-22

- 25 lbs/11.3 kg counterweight kit

Add-on Counterweight Kit

RST-0007-29

- 10 lb/4.5 kg add on counterweight kit for VHRC Series



3875 Cypress Drive
Petaluma, California 94954-5635
United States
800-228-2555
707-773-1100 outside US
707-773-1180 fax

Postbus 10462
5000 JL Tilburg
Netherlands
+31 (0) 88 627 2600
+31 (0) 88 627 2609 fax

No. 7, Lane 116, Wugong 3rd Road
Wugu Township, Taipei County 248
Taiwan (R.O.C.)
+01 707-773-1100
+01 707-773-1180 fax

gcx.com
one tool. one source.