

# GCX<sup>®</sup> Tablet Roll Stand

## Today's solution for tomorrow's telehealth needs

- Equipped with a fixed or variable height tablet arm for optimal viewing and access for video conferencing, electronic medical records, and patient services from anywhere within a hospital setting.
- From translation apps and online patient registrations to virtual rounding and telemedicine, our GCX Tablet Roll Stand allows time-pressed nurses, clinicians, and physicians to meet with patients and specialists no matter their location.
- Optional accessories such as speaker mounts and bins for additional storage.


*Note: Tablet enclosure sold separately. Tablet not provided.*



### FEATURES AND BENEFITS



**Optional PRO-ADJUST™ tablet arm** allows the user to reposition the tablet while remaining firmly secured to the roll stand, reducing the risks of accidental drops or removals.



**Power box** to secure an AC adapter or USB hub along the channel to prevent accidental or intentional unplugging.



**Stable, 18"-by-18"** base and four-inch locking casters seamlessly move the roll stand in and out of a room or viewing area with ease.



**Upto 100% of the cables** are concealed for better infection control.

### ON-SITE SERVICES

- GCX On-Site Services (OSS) is a seamless, cost-effective and expert service for deploying and maintaining GCX products.
- Click [here](#) for more information on GCX On-Site Services.

**Are you ready to roll? Then contact us to learn how the GCX Table Roll Stand can improve your provider and patient experiences to remaining competitive.**

[www.gcx.com](http://www.gcx.com) | 800-228-2555

GCX quality management processes are ISO 13485:2003 certified. GCX, PolyMount, and PolyQuip are registered trademarks, and VHRC, VHC, VHM, VHRS, PRC, M Series, and Stor-Lock are trademarks, of GCX Corporation. All other brand names, product names, or trademarks belong to their respective holders. © Copyright 2021.

Updated 01-11-2024

