


VHM-25™ Counter Top Mount

VHM-25 Counter Top Mount

- 4" x 4"/10.2 cm x 10.2 cm base with 3"/7.6 cm square hole pattern
- Small profile - space saving design
- 7 - 20 lb / 3.2 - 9.1 kg Weight limit
- Tilt/swivel adjustment
- Portrait/landscape adjustment
- 100/75 mm VESA compatible
- Smooth continuous surface with minimal seams for easy cleaning
- Height adjustability 12"/30.5 cm
- 7"/17.8 cm Angled extension for optimal reach
- Flush fold-over ability for minimal foot print
- 20"/50.8 cm maximum reach
- Product weight: 7 lbs/3.2 kg

Note: Patent Pending

WS-0008-05



3875 Cypress Drive
Petaluma, California 94954-5635
United States
800-228-2555
707-773-1100 outside US
707-773-1180 fax

Postbus 10462
5000 JL Tilburg
Netherlands
+31 (0) 88 627 2600
+31 (0) 88 627 2609 fax

No. 7, Lane 116, Wugong 3rd Road
Wugu Township, Taipei County 248
Taiwan (R.O.C.)
+01 707-773-1100
+01 707-773-1180 fax

gcx.com
one tool. one source.