

# VHM-P

VARIABLE HEIGHT ARM

The industry's leading adjustable mounting arm delivers unparalleled usability.

PATENTED



## FEATURES

### Cable Management

Enables easier cleaning, provides simple access for service, and creates a professional-looking installation.

- Up to eight (8) cables can be accommodated
- 80 percent (approx.) of cabling is concealed
- Covers can be easily removed and reinstalled

### Precision

The arm moves smoothly into place at every point across its full range of motion for optimal positioning and viewing.

- Seamless, unibody construction
- Swivel/tilt and pivot-tension adjustments and parallel linkage for consistent viewing angle

### Durability and Cleanability

Withstands constant use and cleaning regimens to outlast generations of devices. Easy to keep clean to fight the spread of infection.

- Medical-grade arm
- Strong, seamless, chemically resistant materials including die-cast aluminum
- Enclosed cables and pivot points



OPTIONAL HANDLE



CABLE MANAGEMENT



SEAMLESS, UNIBODY CONSTRUCTION



ENCLOSED CABLES AND PIVOT POINTS

## TECHNICAL SPECIFICATIONS

- **Vertical range:** 16" / 40.6 cm total travel
- **Lateral range:** +/- 90 degrees
- **Load range:** VHM-P: 15 - 30 lbs / 6.8 - 13.6 kg or 25 - 45 lbs / 11.3 - 20.4 kg - Max Mounting Plate Load: 30 lbs / 13.6 kg
- **Product weight:** 9.4 lb / 4.3 kg
- Available with 8" / 20.3 cm or 14" / 35.6 cm Extension

## AVAILABLE ACCESSORIES

- Cable Management Hooks
- Down Post for Ancillary Devices
- Module Housing Hanger
- Non-Locking Handle

Channel sold separately.

