



FOR IMMEDIATE RELEASE

Smaller VHM-25™ support arm heralds new era in design for GCX mounting solutions

PETALUMA, Calif. (Feb. 13, 2012) – GCX Corporation, the worldwide leader in medical instrument and IT mounting solutions, introduces the new VHM-25™ variable height support arm, a small-profile version of its best-selling VHM™ Series of Variable Height Mounts.

The new VHM-25 is available in two models: a channel mount version and a countertop “bolt-down” version with hidden fasteners. Horizontal and angled 7-inch (17.8-cm) extension options provide longer reach, and hide cables internally for improved cleanability.

GCX created the space-saving design specifically for health care environments. The small amount of space the VHM-25 occupies makes it a natural solution space planners can use to easily integrate flat-screen displays into tight settings.

The VHM-25 is not only smaller, but more adjustable. It gives the mounted device 12 inches (30.5 cm) of height adjustment, side-to-side pivot, tilt angle adjustment, and rotation between portrait and landscape orientations.

The next-generation design of the VHM-25 offers health care institutions a support arm that’s both easy to clean and durable. The seamless top surface is impervious to liquids, while the cast-aluminum housing will stand up to rough handling and cleaning solutions normally used in health care facilities.

“The new VHM-25 support arm is more than an aesthetic update for GCX”, said Vice President of Development Cris Daugbjerg.

“This latest development from GCX points the direction for next-generation support arms that must meet health care institutions’ long-term mounting needs,” he said. “Our customers expect GCX products to embody solid construction and an upgradeable, open-architecture design approach. The VHM-25 now introduces these key GCX attributes to a small form factor.”

###

Note to Editors: High-resolution images of the VHM-25 are available by clicking “Download” on your selection(s) under “View More Images” at the following links.

Channel Mount

<http://www.gcx.com/products/vhm-25-channel-mount>

Counter Top Mount

<http://www.gcx.com/products/vhm-25-counter-top-mount>

About GCX Corporation

GCX Corporation has been serving the health care industry since 1971, providing medical instrument and IT mounting solutions. Products are developed exclusively for the health care field and are marketed directly to hospitals. GCX also has custom product development relationships with Original Equipment Manufacturers. Major product lines include wall mounts, roll stands, ceiling mounts, countertop mounts, pole mounts and more, along with a variety of mounting accessories. GCX mounting solutions feature space-saving, ergonomic designs to improve both equipment and patient access.

More Information

For more information, please contact Kevin Merritt, Director of Marketing, GCX Corporation; 707.773.1100; kevin@gcx.com